


How To: Install g.MOBIIab+ Driver Version 1.14.01

Perform the following steps for installation

1. If there is any old version of the g.MOBIIab+ driver package on your computer please uninstall it.
2. It is highly recommended to turn off the User Account Control (UAC) of the Windows 7 operating system. Please see last page of this document how to do this.
3. Close all running applications.
4. Insert the g.tec installation CD, open the g.MOBIIab/g.MOBIIab+ Driver directory, select the correct directory for the architecture of the PC (Win32 or Win64). To install the g.MOBIIab+ Driver double-click the `setup.exe`. This opens the following welcome dialog (see picture below):



5. Click **Next**. Please read the **License Agreement** for g.MOBIIab+ Driver and if you agree with the terms click **I Agree** and **Next**. Then just follow the steps on the screen.

Turn off User Account Control (UAC)

Please switch off the User Account Control Settings on your system: Click on the Windows **Start** button → **Control Panel** → **System and Security** → **Change User Account Control Settings**. Set the slider to the lowest position (Never notify).

