g.tec medical engineering GmbH

updated: July. 3rd, 2013

## **Remote Dongle Update**

1) Save the file called "gtecHSAPHLupdate.exe" included in this email or

2) Transfer the "gtecHSAPHLupdate.exe" to the PC, where the dongle is used or at least the hasp-hl drivers are installed

v1.12.00

- 3) Plug in the green hasp-hl dongle, that should be up graded check the red LED shining, showing correct function of the dongle
- 4) Remove all other dongles from the PC
- 5) Start "gtecHASPHLupdate.exe" -> select the tab called "Collect Key Status Information"
- 6) Press the "Collect information" button and select where to save the collected info wait some seconds
- 7) Email the saved file (\*.c2v) to g.tec (support@gtec.at) or to your contact person at g.tec.
- 8) We will prepare the upgrade and you will get back the updated file (\*.v2c) via email.
- 9) Start again the program and select the tab called "Apply License Update"
- 10) Select the file you got from us (\*.v2c) and press the "Apply update" button
- 11) Your dongle will be updated